

HIF1B/Arnt Rabbit Polyclonal Antibody For If/Icc, Elisa - P27540

Item Number: CM0024193



Introduction

HIF1 β /ARNT rabbit polyclonal antibody (pAb) validated for IF/ICC and ELISA. Detects human, mouse, rat ARNT (SWISS: P27540). Ideal for hypoxia response and aryl hydrocarbon receptor signaling research.

[Learn More](#)

Parameter	Value
Product Name	HIF1 β /ARNT Rabbit pAb
Remarks/Alias	HIF1B; TANGO; bHLHe2; HIF1BETA; HIF-1beta; HIF1-beta; HIF-1-beta; HIF1 β /ARNT
Species	Human
GeneID (Human)	405
GeneID	405
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 465-570 of human HIF1 β /ARNT (NP_001659.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27540
Protein Weight	87kDa
Shipping	Ice bag