

High-Performance Rabbit Monoclonal Antibody Panel For Human Treg Flow Cytometry - Human Treg Panel

Item Number: CM0011940

Introduction

A rabbit monoclonal antibody panel for human regulatory T cell (Treg) analysis by flow cytometry. Enables precise immunophenotyping of Treg populations, supporting immunology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Human Treg Panel
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
Shipping	Ice bag