

Histone H1T Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q07133

Item Number: CM0023318

Introduction

High-quality rabbit polyclonal antibody against testis-specific Histone H1t, validated for WB, IF-P, and ELISA in mouse and rat samples—ideal for chromatin and meiosis IVD research applications.

[Learn More](#)

Parameter	Value
Product Name	Histone H1t Rabbit pAb
Remarks/Alias	H1t; H1-6; H1.6; H1ft; Hist1h1t; Histone H1t
Species	Mouse
GeneID (human)	H1
GeneID	107970
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q07133
Protein Weight	22kDa
Shipping	Ice bag