

## Hla-A Rabbit Mab - P01891

Item Number: CM0002260

### Introduction

Rabbit monoclonal antibody targeting human HLA-A, validated for WB, IHC-P, ELISA. Cross-reacts with human samples, predicted molecular weight 41 kDa. Suitable for immunology and transplantation research.

[Learn More](#)

| Parameter        | Value                                                                                            |
|------------------|--------------------------------------------------------------------------------------------------|
| Product Name     | HLA-A Rabbit mAb                                                                                 |
| Remarks/Alias    | HLAA; HLA-A                                                                                      |
| Species          | Human                                                                                            |
| Gene ID          | 3105                                                                                             |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HLA-A (P04439) |
| Source           | Rabbit                                                                                           |
| Category         | Monoclonal Antibodies                                                                            |
| Application      | WB, IHC-P, ELISA                                                                                 |
| Cross Reactivity | Human                                                                                            |
| SWISS            | P01891                                                                                           |
| Protein Weight   | 41kDa                                                                                            |
| Shipping         | Ice bag                                                                                          |