

Hla-Dpb1 Rabbit Pab For Wb, Ihc-P, Elisa - P04440

Item Number: CM0013130



Introduction

Rabbit polyclonal antibody against human HLA-DPB1, a MHC class II antigen presenting exogenous antigens to CD4 T cells. Validated for WB, IHC-P, ELISA with cross-reactivity in human, mouse, rat.

[Learn More](#)

Specification	Details
Product Name	HLA-DPB1 Rabbit pAb
Remarks/Alias	DPB1; HLA-DP; HLA-DPB; HLA-DP1B; HLA-DPB1
Species	Human
Gene ID	3115
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-225 of human HLA-DPB1 (NP_002112.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04440
Protein Weight	29 kDa
Shipping	Ice bag