

Hla-Drb1 Rabbit Polyclonal Antibody For Wb And Elisa - P01911 / P01912 / Q29974

Item Number: CM0030315



Introduction

High-quality HLA-DRB1 rabbit polyclonal antibody for Western blot and ELISA applications. Detects human and mouse MHC class II beta chain (30 kDa). Ideal for immunology research on antigen presentation and T cell responses.

[Learn More](#)

Parameter	Specification
Product Name	HLA-DRB1 Rabbit pAb
Remarks/Alias	SS1; DRB1; HLA-DRB; HLA-DR1B; HLA-DRB1
Species	Human
GeneID (Human)	3123
GeneID	3123
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 30-230 of human HLA-DRB1 (NP_002115.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P01911 / P01912 / Q29974
Protein Weight	30kDa
Shipping	Ice bag