

Hla-G Rabbit Polyclonal Antibody For Wb, Elisa - P17693

Item Number: CM0029110



Introduction

Rabbit polyclonal antibody against human HLA-G, validated for Western blot and ELISA. Cross-reacts with mouse and rat. Useful for studying immune tolerance at the maternal-fetal interface, reproductive immunology, and immune checkpoint research. Predicted molecular weight 38 kDa.

[Learn More](#)

Parameter	Value
Product Name	HLA-G Rabbit pAb
Remarks/Alias	MHC-G; HLA-G
Species	Human
GeneID (Human)	3135
GeneID	3135
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 20-230 of human HLA-G (NP_002118.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17693
Protein Weight	38kDa
Shipping	Ice bag