

Hnrnpf Rabbit Polyclonal Antibody For Wb, Elisa - P52597

Item Number: CM0023337



Introduction

HNRNPF rabbit polyclonal antibody validated for WB and ELISA; reacts with human and mouse. Suitable for studying alternative splicing and RNA-binding protein function.

[Learn More](#)

Parameter	Value
Product Name	HNRNPF Rabbit pAb
Remarks/Alias	HNRNPF; mcs94-1; OK/SW-cl.23; HNRNPF
Species	Human
Gene ID (Human)	3185
Gene ID	3185
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 315-415 of human HNRNPF (NP_004957.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P52597
Protein Weight	46kDa
Shipping	Ice bag