

Hnrnpm Rabbit Mab For Wb, Ip, Ihc-P, If/Icc, Elisa - P52272

Item Number: CM0003367



Introduction

HNRNPM Rabbit mAb is a rabbit monoclonal antibody recognizing HNRNPM (78 kDa) in human, mouse, and rat. Suitable for WB, IP, IHC-P, IF/ICC, and ELISA, it supports research into pre-mRNA splicing and carcinoembryonic antigen signaling.

[Learn More](#)

Parameter	Value
Product Name	HNRNPM Rabbit mAb
Remarks/Alias	CEAR; HNRPM; HTGR1; NAGR1; HNRPM4; HNRNPM4; hnRNP M; HNRNPM
Species	Human
GeneID (Human)	4670
GeneID	4670
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human HNRNPM (P52272)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52272
Protein Weight	78 kDa
Shipping	Ice bag