


Hnrnpr Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O43390

Item Number: CM0020040



Introduction

HNRNPR rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat heterogeneous nuclear ribonucleoprotein R, a key player in mRNA processing. Ideal for RNA biology research.

[Learn More](#)

Parameter	Value
Product Name	HNRNPR Rabbit pAb
Remarks/Alias	HNRPR; NEDDFS; hnRNP-R; HNRNPR
Species	Human
GeneID (Human)	10236
GeneID	10236
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 550-633 of human HNRNPR (NP_005817.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43390
Protein Weight	71kDa
Shipping	Ice bag