

Homer1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q86Ym7

Item Number: CM0028642

Introduction

Rabbit polyclonal antibody to HOMER1. Validated for Western blot, IF-P, and ELISA. Cross-reacts with human, mouse, and rat. Scaffold protein involved in synaptic signaling and plasticity. Ideal for neuroscience research on metabotropic glutamate receptor coupling and dendritic spine morphology.

[Learn More](#)

Parameter	Value
Product Name	HOMER1 Rabbit pAb
Remarks/Alias	HOMER; SYN47; Ves-1; HOMER1A; HOMER1B; HOMER1C; HOMER1
Species	Human
GeneID (Human)	9456
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 125-354 of human HOMER1 (NP_004263.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86YM7
Protein Weight	40kDa
Shipping	Ice bag