

Hopx Rabbit Polymab® For Wb, Ip, If-P, Ihc-P, Elisa - Q9Bpy8

Item Number: CM0002214

Introduction

Recombinant monoclonal antibody targeting human HOPX (Homeodomain-only protein). Validated for WB, IP, IF-P, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for cardiac development and cancer research.

[Learn More](#)

Parameter	Value
Product Name	HOPX Rabbit PolymAb®
Remarks/Alias	HOD; HOP; OB1; LAGY; TOTO; CAMEO; NECC1; SMAP31
Species	Human
GeneID (Human)	84525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-73 of human HOPX (NP_631957.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BPY8
Protein Weight	8 kDa
Shipping	Ice bag