

Hoxb4 Rabbit Monoclonal Antibody For Wb And Elisa - P17483

Item Number: CM0002241

Introduction

Recombinant rabbit monoclonal antibody against human HOXB4, validated for Western blot and ELISA. Suitable for studying homeobox transcription factors in developmental biology and gene regulation research.

[Learn More](#)

Parameter	Value
Product Name	HOXB4 Rabbit mAb
Remarks/Alias	HOX2; HOX2F; HOX-2.6; HOXB4
Species	Human
Gene ID	3214
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human HOXB4 (P17483)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P17483
Protein Weight	28kDa
Shipping	Ice bag