

Hoxb5 Rabbit Pab For Wb, If/Icc, Elisa - P09067

Item Number: CM0027933

Introduction

High-quality HOXB5 rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat species. Ideal for developmental biology and transcription research. Approximately 29 kDa.

[Learn More](#)

Parameter	Value
Product Name	HOXB5 Rabbit pAb
Remarks/Alias	HOX2; HU-1; HOX2A; Hox2.1; HHO.C10; HOXB5
Species	Human
Gene ID (Human)	3215
Gene ID	3215
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 170-269 of human HOXB5 (NP_002138.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09067
Protein Weight	29kDa
Shipping	Ice bag