

Hrsp12 Rabbit Polyclonal Antibody For Wb, Elisa - P52758

Item Number: CM0026424



Introduction

High-quality rabbit polyclonal antibody targeting human HRSP12 (RIDA), validated for WB and ELISA. Suitable for studying enamine/imine metabolism and m6A-dependent RNA cleavage. Cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	HRSP12 Rabbit pAb
Remarks/Alias	PSP; P14.5; UK114; HRSP12; hp14.5
Species	Human
Gene ID	10247
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-137 of human HRSP12 (NP_005827.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52758
Protein Weight	14kDa
Shipping	Ice bag