

Hs2St1 Rabbit Pab - Q7Lga3

Item Number: CM0015117

Introduction

HS2ST1 Rabbit pAb is a polyclonal antibody validated for WB and ELISA. It detects human, mouse, and rat HS2ST1, a 42kDa heparan sulfate 2-O-sulfotransferase involved in kidney development and HS biosynthesis.

[Learn More](#)

Parameter	Value
Product Name	HS2ST1 Rabbit pAb
Remarks/Alias	NFSRA; dj604K5.2; HS2ST1
Species	Human
GeneID (Human)	9653
GeneID	9653
Immunogen	Recombinant protein/Recombinant fusion protein containing a sequence corresponding to amino acids 1-229 of human HS2ST1 (NP_001127964.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7LGA3
Protein Weight	42kDa
Shipping	Ice bag