

Hsc70/Hspa8 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - P11142

Item Number: CM0020556



Introduction

Rabbit polyclonal antibody against human Hsc70/HSPA8 (Heat shock cognate 71 kDa protein, HSPA8). Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse, rat, and monkey. Suitable for chaperone, autophagy, and stress research.

[Learn More](#)

Parameter	Value
Product Name	Hsc70/HSPA8 Rabbit pAb
Remarks/Alias	LAP1; HSC54; HSC70; HSC71; HSP71; HSP73; LAP-1; NIP71; HEL-33; HSPA10; HEL-S-72p; Hsc70/HSPA8
Species	Human
GeneID (Human)	3312
GeneID	3312
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 510-615 of human Hsc70/HSPA8 (NP_006588.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	P11142
Protein Weight	71kDa
Shipping	Ice bag