


Hsp27/Hspb1 Rabbit Pab - P04792

Item Number: CM0012905



Introduction

HSP27/HSPB1 Rabbit Polyclonal Antibody validated for Western blot (WB), Immunohistochemistry (IHC-P), and ELISA. Reacts with human, mouse, and rat HSP27/HSPB1 (SWISS P04792). Targets heat shock protein beta-1, a molecular chaperone involved in stress response.

[Learn More](#)

Parameter	Value
Product Name	HSP27/HSPB1 Rabbit pAb
Remarks/Alias	CMT2F; HMN2B; HSP27; HSP28; Hsp25; SRP27; HS.76067; HEL-S-102; HSP27/HSPB1
Species	Human
GeneID (Human)	3315
GeneID	3315
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human HSP27/HSPB1 (NP_001531.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04792
Protein Weight	23kDa
Shipping	Ice bag