

Hspa14 Rabbit Mab For Wb, Ihc-P, If/Icc, Elisa - Q0Vdf9

Item Number: CM0006620

Introduction

Recombinant rabbit monoclonal antibody against human HSPA14 (HSP70-4, HSP70L1) for WB, IHC-P, IF/ICC, ELISA. Validated in Human and Mouse. Cytosolic chaperone involved in nascent polypeptide folding.

[Learn More](#)

Parameter	Value
Product Name	HSPA14 Rabbit mAb
Remarks/Alias	HSP70-4; HSP70L1; MSANTD7; HSPA14
Species	Human
Gene ID (Human)	51182
Gene ID	51182
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HSPA14 (Q0VDF9).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q0VDF9
Protein Weight	55kDa
Shipping	Ice bag