


Hspb11 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y547

Item Number: CM0023722



Introduction

Rabbit polyclonal antibody against HSPB11 (IFT25, PP25) validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for ciliary biology, SHH signaling, and fertility research.

[Learn More](#)

Parameter	Value
Product Name	HSPB11 Rabbit pAb
Remarks/Alias	PP25, FAP232, HSPB11, C1orf41, CFAP232, HSPCO34
Species	Human
Gene ID (Human)	51668
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-144 of human HSPB11 (NP_057210.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y547
Protein Weight	16kDa
Shipping	Ice bag