


Htatsf1 Rabbit Pab - O43719

Item Number: CM0029709



Introduction

Rabbit polyclonal antibody against HTATSF1, a spliceosome component and HIV Tat cofactor. Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for nuclear splicing and transcription research.

[Learn More](#)

Parameter	Value
Product Name	HTATSF1 Rabbit pAb
Remarks/Alias	TATSF1; TAT-SF1; dj196E23.2; HTATSF1
Species	Human
GeneID (Human)	27336
GeneID	27336
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human HTATSF1 (NP_055315.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43719
Protein Weight	86kDa
Shipping	Ice bag