


Htra3 Rabbit Polyclonal Antibody For Wb And Elisa - P83110

Item Number: CM0013932



Introduction

Rabbit polyclonal antibody targeting human HTRA3 (serine protease, tumor suppressor, placental development). Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for reproductive biology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	HTRA3 Rabbit pAb
Remarks/Alias	Prsp; Tasp; HTRA3
Species	Human
Gene ID (Human)	94031
Gene ID	94031
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 334-453 of human HTRA3 (NP_444272.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P83110
Protein Weight	49kDa
Shipping	Ice bag