

Human Activated T Cell Panel - Rabbit Monoclonal Antibody

Item Number: CM0009186

Introduction

Human Activated T Cell Panel, rabbit monoclonal antibodies for flow cytometry analysis of human activated T cells. Designed for research use in immunology and cell biology.

[Learn More](#)

Parameter	Details
Product Name	Human Activated T Cell Panel
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Human
Shipping	Ice bag