

Human Ager Rabbit Pab For Wb, Elisa - Q15109

Item Number: CM0019034



Introduction

Polyclonal antibody targeting human AGER (RAGE) validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Suitable for inflammation, diabetes, and neurodegenerative disease research.

[Learn More](#)

Parameter	Value
Product Name	Human AGER Rabbit pAb
Remarks/Alias	RAGE; sRAGE; SCARJ1; Human AGER
Species	Human
Gene ID (Human)	177
Gene ID	177
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-344 of human AGER (NP_001127.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q15109
Protein Weight	43kDa
Shipping	Ice bag