

Human Ascs B Cell Panel - Human Ascs B Cell Panel

Item Number: CM0010954

Introduction

Human ASCs B cell monoclonal antibody panel designed for flow cytometry (FC) applications, specifically reactive with human samples. Ideal for research into antibody-secreting B cell identification and characterization.

[Learn More](#)

Parameter	Value
Product Name	Human ASCs B cell Panel
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Human
Shipping	Ice bag