


Human IgA1 Rabbit Monoclonal Antibody For Wb And Elisa - P01876

Item Number: CM0004054



Introduction

CamelBio offers this high-quality rabbit monoclonal antibody targeting human IgA1, validated for WB and ELISA applications. Suitable for research on mucosal immunity, immunoglobulin function, and B lymphocyte biology.

[Learn More](#)

Parameter	Value
Product Name	Human IgA1 Rabbit mAb
Remarks/Alias	IgA1
Species	Human
Gene ID	3493
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 99-353 of human IgA1 (P01876)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P01876
Protein Weight	37KDa/42KDa
Shipping	Ice bag