

Human IgE Rabbit Mab For Wb, Fc, Elisa - P01854

Item Number: CM0006982



Introduction

Monoclonal rabbit antibody targeting human IgE heavy chain constant region (IGHE, P01854). Validated in WB, FC, ELISA. Recognizes secreted and membrane IgE, essential for allergy, inflammation, and B cell research.

[Learn More](#)

Parameter	Description
Product Name	Human IgE Rabbit mAb
Remarks/Alias	Igh2; IGHE; IGHEP1; Human IgE
Species	Human
Gene ID	3497 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 104-428 of human IgE (P01854).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	P01854
Protein Weight	47kDa
Shipping	Ice bag