


Human IgG (Fc) Rabbit Monoclonal Antibody For Wb, Ihc-P, If-P, Elisa - P01857 / P01859 / P01860 / P01861

Item Number: CM0006024



Introduction

High-quality Human IgG (Fc) Rabbit Monoclonal Antibody, validated for WB, IHC-P, IF-P, ELISA. Recognizes human IgG subclasses (P01857, P01859, P01860, P01861). Recombinant immunogen, 52kDa. Ideal for immunological and diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	Human IgG (Fc) Rabbit mAb
Species	Human
GeneID (Human)	3,500,350,135,023,500
GeneID	3,500,350,135,023,500
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-330 of human IgG (Fc) mAb (P01857).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P01857, P01859, P01860, P01861
Protein Weight	52kDa
Shipping	Ice bag