


Human IgG1/Immunoglobulin Heavy Constant Gamma 1 Rabbit Polyclonal Antibody For Wb, Elisa - P01857

Item Number: CM0024143



Introduction

Rabbit polyclonal antibody targeting human IgG1 (heavy constant gamma 1, IGHG1). Suitable for Western blotting (WB) and ELISA. Cross-reacts with human and mouse. Immunogen corresponds to amino acids 21-130 (CH1 domain).

[Learn More](#)

Parameter	Value
Product Name	Human IgG1/Immunoglobulin heavy constant gamma 1 Rabbit pAb
Species	Human
GeneID (Human)	3500
GeneID	3500
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-130 of human IgG1/Immunoglobulin heavy constant gamma 1 (P01857)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P01857
Protein Weight	44kDa
Shipping	Ice bag