


## Human IgM Rabbit Pab For Wb, Elisa - P01871

Item Number: CM0026493



### Introduction

Human IgM Rabbit pAb targets the human IgM heavy chain constant region (UniProt P01871). Validated for Western blot and ELISA, this polyclonal antibody is essential for research on humoral immunity, B cell biology, and immunoglobulin structure-function studies.

[Learn More](#)

| Parameter        | Value                                                                                                        |
|------------------|--------------------------------------------------------------------------------------------------------------|
| Product Name     | Human IgM Rabbit pAb                                                                                         |
| Remarks/Alias    | MU; VH; AGM1; Human IgM                                                                                      |
| Species          | Human                                                                                                        |
| GeneID (Human)   | 3507                                                                                                         |
| GeneID           | 3507                                                                                                         |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 213-453 of human IGHM (P01871) |
| Source           | Rabbit                                                                                                       |
| Category         | Polyclonal Antibodies                                                                                        |
| Application      | WB, ELISA                                                                                                    |
| Cross Reactivity | Human                                                                                                        |
| SWISS            | P01871                                                                                                       |
| Protein Weight   | 49kDa                                                                                                        |
| Shipping         | Ice bag                                                                                                      |