

Human Neutrophils Panel Monoclonal Antibody For Flow Cytometry - Human Neutrophils Panel

Item Number: CM0008196

Introduction

High-quality monoclonal antibody panel targeting human neutrophils, validated for flow cytometry (FC). Ideal for immunophenotyping and research on neutrophil biology. Supplied with ice bag shipping.

[Learn More](#)

Parameter	Value
Product Name	Human Neutrophils Panel
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
Shipping	Ice bag