

Human T Cells Panel - Rabbit Monoclonal Antibody

Item Number: CM0012151

Introduction

Rabbit monoclonal antibody panel for human T cell analysis by flow cytometry (FC). Designed for immuno-oncology and T cell research applications. Cross-reacts with human samples. Shipped on ice bag.

[Learn More](#)

Parameter	Value
Product Name	Human T Cells Panel
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
Shipping	Ice bag