

Human Tbnk Panel Monoclonal Antibody Cocktail - Human Tbnk Panel

Item Number: CM0011898

Introduction

Human TBNK Panel consisting of rabbit monoclonal antibodies for flow cytometry detection of T, B, and NK cells in human samples. Suitable for immunophenotyping in research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Human TBNK Panel
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
Shipping	Ice bag