

Human Thioredoxin Reductase 1/Trxr1 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - Q16881

Item Number: CM0021334



Introduction

Rabbit polyclonal antibody against human Thioredoxin Reductase 1 (TXNRD1) for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Suitable for redox biology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Human Thioredoxin Reductase 1/TRXR1 Rabbit pAb
Remarks/Alias	TR; TR1; TXNR; TRXR1; GRIM-12; Thioredoxin reductase 1 (TXNRD1)
Species	Human
Gene ID	7296
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 410-490 of human Thioredoxin reductase 1 (TXNRD1) (NP_001087240.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16881
Protein Weight	71kDa
Shipping	Ice bag