

Human Tlr3 Rabbit Pab - O15455

Item Number: CM0022877



Introduction

High-quality Human TLR3 Rabbit pAb for WB, IF/ICC, ELISA. Recognizes human TLR3, a key innate immune receptor that detects viral dsRNA and activates inflammatory responses. Suitable for immunology and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	Human TLR3 Rabbit pAb
Remarks/Alias	TLR3; CD283; IIAE2; toll-like receptor 3; Human TLR3
Species	Human
Gene ID (Human)	7098
Gene ID	7098
Immunogen	Recombinant protein corresponding to a sequence within amino acids 24-339 of human TLR3 (NP_003256.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O15455
Protein Weight	72kDa/103kDa
Shipping	Ice bag