


Hydroxyl-Histone H2A-Y39 (P04908/P0C0S8) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa

Item Number: CM0000072



Introduction

Rabbit monoclonal antibody specific for hydroxylated histone H2A at tyrosine 39 (P04908/P0C0S8), validated for Western blot, IHC-P, and ELISA. Reacts with human, mouse, rat, and other species. Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Hydroxyl-Histone H2A-Y39 Rabbit mAb
Remarks/Alias	H2A.1; H2A.2; H2A/a; H2AC4; H2AFA; HIST1H2AE; Hydroxyl-Histone H2A-Y39
Species	Human
GeneID (human)	3012
GeneID	3012/8329
Immunogen	Synthetic Peptide: A synthetic hydroxylated peptide around Y39 of human Histone H2A (P04908)
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P04908/P0C0S8
Protein Weight	14kDa
Shipping	Ice bag