

Icam-1/Cd54 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P13597

Item Number: CM0005476

Introduction

ICAM-1/CD54 Rabbit Monoclonal Antibody validated for WB, IF-P, IHC-P, ELISA. Reacts with mouse and rat. Essential for studying leukocyte adhesion and transmigration.

[Learn More](#)

Parameter	Value
Product Name	ICAM-1/CD54 Rabbit mAb
Remarks/Alias	CD54; Ly-47; Icam-1; MALA-2; ICAM-1/CD54
Species	Mouse
GeneID (Human)	3383
GeneID	15894
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-485 of mouse ICAM-1/CD54 (NP_034623.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P13597
Protein Weight	59kDa
Shipping	Ice bag