

Id4 Rabbit Polyclonal Antibody For Wb, Elisa - P47928

Item Number: CM0018843

Introduction

Rabbit polyclonal antibody targeting human ID4, a transcriptional regulator implicated in cell growth, differentiation, and cancer. Validated for WB and ELISA, with cross-reactivity to mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	ID4 Rabbit pAb
Remarks/Alias	IDB4; bHLHb27; ID4
Species	Human
GeneID	3400 (Human)
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human ID4 (NP_001537.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P47928
Protein Weight	17kDa
Shipping	Ice bag