

Ifluor 594/Abflo 594 Rabbit Monoclonal Antibody For Db And Elisa - If594

Item Number: CM0010413

Introduction

Ifluor 594/ABflo 594 Rabbit Monoclonal Antibody (iF594) is a species-independent antibody for DB and ELISA applications. It specifically recognizes the Ifluor 594 dye, making it an essential tool for fluorescent dye detection in immunoassays.

[Learn More](#)

Parameter	Specification
Product Name	Ifluor 594/ABflo 594 Rabbit mAb
Remarks/Alias	iF594
Species	Independent species
Immunogen	Chemical compounds corresponding to Ifluor 594
Source	Rabbit
Category	Monoclonal Antibodies
Application	DB, ELISA
Cross Reactivity	Independent species
Shipping	Ice bag