


Ifn-Beta Rabbit Monoclonal Antibody For Wb, Elisa - P01575

Item Number: CM0007131



Introduction

Rabbit monoclonal antibody against mouse IFN-beta, suitable for Western blot (WB) and ELISA applications. IFN-beta is a key innate immune cytokine with antiviral and immunoregulatory functions. Ideal for immunology and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	IFN-beta Rabbit mAb
Remarks/Alias	Ifb; IFNB; If1da1; IFN-beta
Species	Mouse
GeneID (Human)	3456
GeneID	15977
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-182 of mouse IFN-beta (NP_034640.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P01575
Protein Weight	22kDa
Shipping	Ice bag