

Ifn Gamma Rabbit Polymab® For Wb And Elisa - P01579 / P01581

Item Number: CM0005598



Introduction

IFN gamma Rabbit PolymAb® (SWISS P01579/P01581), a monoclonal antibody detecting human and rat interferon gamma via WB and ELISA. Ideal for immunological and infection research studies.

[Learn More](#)

Parameter	Value
Product Name	IFN gamma Rabbit PolymAb®
Remarks/Alias	IF2f; IFNG2
Species	Human/Rat
GeneID (Human)	3458
GeneID	3458/25712
Immunogen	Recombinant protein of human/rat IFN gamma
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P01579/P01581
Protein Weight	18kDa
Shipping	Ice bag