

Ifnar2 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - P48551

Item Number: CM0019734



Introduction

Rabbit polyclonal antibody targeting human IFNAR2 (IFNARB, IFNABR) validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying type I interferon signaling and JAK-STAT pathway activation.

[Learn More](#)

Parameter	Value
Product Name	IFNAR2 Rabbit pAb
Remarks/Alias	IFN-R; IMD45; IFNABR; IFNARB; IFN-R-2; IFN-alpha-REC; IFNAR2
Species	Human
GeneID (Human)	3455
GeneID	3455
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-243 of human IFNAR2 (NP_997468.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48551
Protein Weight	58kDa
Shipping	Ice bag