

Ifnar2 Rabbit Polyclonal Antibody (Ko-Validated) - [Ko Validated] Ifnar2 Rabbit Pab

Item Number: CM0023765

Introduction

Rabbit polyclonal antibody targeting human IFNAR2 (IFN-R, IFNABR, IFNARB). Validated in Knockout (KO) models. Suitable for Western Blot (WB) and ELISA applications. Cross-reacts with human, mouse, and rat. Ideal for interferon signaling and immunology research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] IFNAR2 Rabbit pAb
Remarks/Alias	IFN-R; IMD45; IFNABR; IFNARB; IFN-R-2; IFN-alpha-REC; R2
Species	Human
Gene ID	3455
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 101-243 of human IFNAR2 (NP_997468.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
Protein Weight	27kDa/37kDa/57kDa
Shipping	Ice bag