

Igf2Bp1/Imp1 Rabbit Mab For Wb, Ihc-P, If/Icc, Ip, Elisa - Q9Nzi8

Item Number: CM0011016



Introduction

Recombinant rabbit monoclonal antibody targeting human IGF2BP1/IMP1 (Q9Nzi8). Validated for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Ideal for studying mRNA transport, translation regulation, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	IGF2BP1/IMP1 Rabbit mAb
Remarks/Alias	IMP1; ZBP1; CRDBP; IMP-1; CRD-BP; VICKZ1
Species	Human
GeneID (Human)	10642
GeneID	10642
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 101-200 of human IGF2BP1/IMP1(NP_006537.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Nzi8
Protein Weight	49kDa/63kDa
Shipping	Ice bag