


Igf2Bp3 Rabbit Polyclonal Antibody (Wb, Ip, If, Ihc-P, Elisa) - O00425

Item Number: CM0026782



Introduction

Rabbit polyclonal antibody targeting human IGF2BP3 (IMP3). Validated for WB, IP, IF/ICC, IHC-P, and ELISA with cross-reactivity to mouse. Useful for studying mRNA transport, stability, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	IGF2BP3 Rabbit pAb
Remarks/Alias	KOC, CT98, IMP3, KOC1, IMP-3, VICKZ3, IGF2BP3
Species	Human
Gene ID	10643
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 300-579 of human IGF2BP3 (NP_006538.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O00425
Protein Weight	64 kDa
Shipping	Ice bag