

Igfbp1 Rabbit Pab For Wb, Elisa - P08833

Item Number: CM0026575



Introduction

Rabbit polyclonal antibody targeting human IGFBP1; validated for WB and ELISA applications. Cross-reacts with mouse and rat. Appropriate for investigating IGF-mediated pathways, cell migration, and apoptotic regulation.

[Learn More](#)

Parameter	Value
Product Name	IGFBP1 Rabbit pAb
Remarks/Alias	AFBP, IBP1, PP12, IGF-BP25, hIGFBP-1, IGFBP1
Species	Human
Gene ID (Human)	3484
Gene ID	3484
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 81-259 of human IGFBP1 (NP_000587.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08833
Protein Weight	28kDa
Shipping	Ice bag