

Igfbp5 Rabbit Pab - P24593

Item Number: CM0016622



Introduction

IGFBP5 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat IGFBP5, a key regulator of IGF signaling and ECM production. Ideal for oncology and bone biology research.

[Learn More](#)

Parameter	Value
Product Name	IGFBP5 Rabbit pAb
Remarks/Alias	IBP5; IGFBP5
Species	Human
Gene ID (Human)	3488
Gene ID	3488
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 68-167 of human IGFBP5 (NP_000590.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24593
Protein Weight	31kDa
Shipping	Ice bag