

Ikb-Z Rabbit Polyclonal Antibody For Ihc-P, Elisa - Q9Byh8

Item Number: CM0021161



Introduction

Rabbit polyclonal antibody targeting human IκB-ζ (NFKBIZ), validated for IHC-P/ELISA with mouse cross-reactivity. Ideal for NF-kappa-B signaling and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	IκB-ζ Rabbit pAb
Remarks/Alias	IKBZ; INAP; MAIL; IκB-zeta; ikappaBzeta; IkappaB-zeta; I-kappa-B-zeta; NFKBIZ
Species	human
GeneID (Human)	64332
GeneID	64332
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human IκB-ζ (NP_113607.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	Q9BYH8
Protein Weight	78kDa
Shipping	Ice bag