

IL-1 Alpha Rabbit Mab - P01582

Item Number: CM0006454

Introduction

Recombinant rabbit monoclonal antibody targeting mouse IL-1 alpha, suitable for Western blot and ELISA. Ideal for studying inflammation, cytokine signaling, and innate immune responses.

[Learn More](#)

Parameter	Value
Product Name	IL-1 alpha Rabbit mAb
Remarks/Alias	IL-1a; IL-1 alpha
Species	Mouse
Gene ID (Human)	3552
Gene ID (Mouse)	16175
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P01582
Protein Weight	31kDa
Shipping	Ice bag