

IL10 Rabbit Pab For Wb And Elisa - P18893

Item Number: CM0020855



Introduction

Mouse IL-10 rabbit polyclonal antibody validated for WB and ELISA. Recognizes mouse IL-10, a key anti-inflammatory cytokine. Recombinant immunogen. Ideal for immune regulation research.

[Learn More](#)

Parameter	Value
Product Name	IL10 Rabbit pAb
Remarks/Alias	CSIF; If2a; IL-10; IL10
Species	Mouse
GeneID (human)	3586
GeneID	16153
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-178 of mouse IL10 (NP_034678.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P18893
Protein Weight	21kDa
Shipping	Ice bag